


STEMI Protocol



ONSET to BALLOON

LOCAL PCI / CATH LABS

Samaritan Hospital (Troy, NY)

- o 518-271-3424 (ED phone)
- o 518-271-3131 (ED fax)
- o Radio 155.400

Albany Medical Center (Albany, NY)

- o 518-262-4660 (Direct phone)
- o 518-262-3194 (Fax)

Dartmouth Hitchcock (Lebanon, NH) (IF USING HELICOPTER ONLY)

- o 877-999-9870 (Direct phone)

SVMC ECG transmission to ED

- o 802-447-5305 (Medical Control)
- o 802-447-2469 (ED fax)
- o EMSECG@phin.org (E-mail)

PEARLS

- Consider ALS intercept
- Consider Aeromedical Transport
 - o Helicopter activation as needed based on clinical scenario and/or cardiogenic shock
- Consider EMS transfer intercept with local paramedic-level squad if current ambulance cannot be taken out of service area for prolonged period of time due to service area coverage needs

TIPS TO PREVENT FALSE ACTIVATION:

- Good skin prep/lead contact
- Correct lead placement
- Level tracing baseline
- Limit motion
- Limit artifact
- Good waveform in all 12 leads